

*Edition 01 – 12/2023*

## *BM-A-5W-14-TX-ENH* *BM-A-5W-14-TX-E-PC* *User Manual*



**Laser Point srl** - Via Burona, 51 - 20055 Vimodrone (Milano) - Italy  
Phone +39 02 27 400 236 - Fax +39 02 25 029 161-[www.LaserPoint.eu](http://www.LaserPoint.eu)


# ***BM-A-5W-14-TX-ENH***

## **User's Manual**

### **Index**

|   |          |
|---|----------|
| <b>1. Safety notes</b>                                    | <b>3</b> |
| <b>2. Introduction to P/N BM-A-5W-14-TX-ENH</b>           | <b>4</b> |
| <b>3. Instrument Specifications</b>                       | <b>4</b> |
| 3.1. Operating Conditions                                 | 4        |
| 3.2. Maximum Ratings                                      | 4        |
| 3.3. Optical and Electrical Specifications                | 5        |
| 3.4. Mechanical Specifications                            | 6        |
| <b>4. Operating Instructions</b>                          | <b>7</b> |
| <b>5. Electrical Connections and Filters</b>              | <b>7</b> |
| <b>6. Warranty</b>  | <b>7</b> |
| <b>7. Limited Liability</b>                               | <b>7</b> |
| <b>8. Compliance to RoHS Directive (RoHS 2015/863)</b>    | <b>8</b> |
| <b>9. European Union WEEE Directive (WEEE 2012/19/EU)</b> | <b>8</b> |

## 1. Safety notes



Before operating this instrument, carefully review the following safety information to avoid personal injury and prevent damage to this instrument or any sensor head connected to it. This instrument does not contain any user-serviceable parts.

The user of Power/Energy measurement instruments must be trained to the use of lasers and their associated risks (ref. EN 60825). LaserPoint is no way liable for any damage resulting from misuse, careless or use above rated limits for the instrument.

The measuring as well as the use of lasers is potentially dangerous. This instrument may operate over wavelengths including non-visible laser radiations.

Proper operating practice in accordance with laser manufacturer's recommendations is crucial; to ensure correct operating procedures, consult the laser manufacturer and your laser safety officer.

Eyewear and other personal protective equipment must be used in compliance with applicable laws and safety regulations.

Be extremely careful with radiation either back-reflected or back-scattered from detector surfaces, housings, mounts and stainless-steel post.

Operate this instrument only within the specified range of operating conditions.

Do not operate this instrument in critical medical environments, in wet or damp conditions or in an explosive atmosphere.

Do not operate this instrument if in suspect of damage or failures. Refer about damaged equipment to Laser Point for qualified service inspection.

## 2. Introduction to P/N BM-A-5W-14-TX-ENH

Blink model BM-A-5W-14-TX-ENH is a Fast response laser sensor for low power lasers, able to measure the energy of single laser pulses with repetition frequency up to 1 MHz for ultra-short pulsed lasers for laser beams with an average power up to 5W.

Its absorber spectral range is broadband, due to the thermal working principle.

## 3. Instrument Specifications

### 3.1. Operating Conditions

| N° | PARAMETER               | CONDITIONS     | SYMB. | MIN | MAX | UNIT |
|----|-------------------------|----------------|-------|-----|-----|------|
| 1  | Environment Temperature |                |       | 10  | 30  | °C   |
| 2  | Relative Humidity       | Non condensing |       | 0   | 90  | %    |

### 3.2. Maximum Ratings

| N° | PARAMETER  | CONDITIONS  | SYMB.    | MAX | UNIT             |
|----|--|---|----------|-----|------------------|
| 1  | Max Power<br>(max 2 sec exposure,<br>max 20% duty cycle) | $E_p=1 \text{ mJ}$ , $\tau_p=200\text{ns}$                    | $P_O$    | 15  | W                |
| 2  | Max Avg Power Density                                    | $E_p=1 \text{ mJ}$ , $\tau_p=200\text{ns}$<br>$P_O=5\text{W}$ | $D_{PO}$ | 0.1 | $\text{kW/cm}^2$ |
| 3  | Max Energy Density                                       | 4ns @1064nm   | $D_t$    | 35  | $\text{mJ/cm}^2$ |
| 4  | Maximum Energy per pulse                                 | 4ns @1064nm,<br>single shot                                   | $D_{Ep}$ | 10  | mJ               |
| 5  | Sensor Max temperature                                   |   |          | 60  | °C               |

### 3.3. Optical and Electrical Specifications


| N° | PARAMETER                                | CONDITIONS  | SYMB.             | MIN          | TYP                               | MAX        | UNIT                           |
|----|--|---|-------------------|--------------|-----------------------------------|------------|--------------------------------|
| 1  | Absorber spectral range                  |   | $\lambda$         | 0.1          | -                                 | 10.6       | $\mu\text{m}$                  |
| 2  | Calibration spectral range*              |   | $\lambda_c$       | -            | 0.355,<br>0.532,<br>1.07;<br>10.6 | -          | $\mu\text{m}$                  |
| 3  | Average Power                            |   | $P_O$             | -            | -                                 | 5          | W                              |
| 4  | Energy Single Shot                       | 1064nm, $\tau_P < 10 \text{ ns}$<br>1064nm, $\tau_P = 10 \mu\text{s}$             | $E_P$             | 0.25<br>20   | -                                 | 40<br>1000 | $\mu\text{J}$<br>$\mu\text{J}$ |
| 5  | Sensitivity                              | $R_L$ (single channel):<br>50 $\Omega$<br>$R_L$ (single channel):<br>1 M $\Omega$ | $S_{1064}$        | -<br>-       | 2<br>4                            | -<br>-     | mV/W<br>mV/W                   |
| 6  | Calibration Accuracy on Energy and Power |   |                   | -            | -                                 | $\pm 5$    | %                              |
| 7  | Measurable Rep. Rate                     | $\tau_P \leq 4 \text{ ns}$  | $f_{\text{rep}}$  | 0**, 1e-3*** | -                                 | 1          | MHz                            |
| 8  | Decay Time                               | $\tau_P \leq 4 \text{ ns}$  | $T_{100-10}$      | -            | 260                               | -          | ns                             |
| 9  | Spatial Uniformity                       | $\phi = 3\text{mm}$ ;<br>80% Active Area  | $\sigma_{xy}$     | -            | -                                 | $\pm 5$    | %                              |
| 10 | Minimum beam diameter                    |   | $\phi$            | 3            | -                                 | -          | mm                             |
| 11 | Power Linearity                          |   | $\sigma_P$        | -            | -                                 | $\pm 2.5$  | %                              |
| 12 | Energy Linearity                         | Vs pulse duration   | $\sigma_{E_P}$    | -            | -                                 | $\pm 2.5$  | %                              |
| 13 | Linearity vs $T_{OP}$                    |   | $\sigma_{T_{OP}}$ | -            | -                                 | $\pm 2.5$  | %                              |
| 14 | Coating Type                             |   | T                 | -            | T                                 | -          | -                              |
| 15 | NEP (with dedicated electronics)         | $B_W = 250 \text{ MHz}$   |                   | -            | 0.1                               | -          | W                              |
| 16 | Diffuser                                 |   | D                 | -            | N                                 | -          | Y/N                            |

\* Other wavelengths on request.

\*\* Bare sensor

\*\*\* With dedicated HSM electronics

### 3.4. Mechanical Specifications



| N° | PARAMETER                | TYP.   | UNIT |
|----|--------------------------|--|------|
| 1  | Outer Dimensions         | 60x60x16.5   | mm   |
| 2  | Weight                   | 130  | g    |
| 3  | Active Area              | 14x14  | mm   |
| 4  | Connection Cable Length  | 2  | m    |
| 5  | Connector Type           | Hirose IX  | -    |
| 8  | Available Absorber Types | Broadband<br>$\lambda = 1064\text{nm} \rightarrow S_{1064}$<br>$\lambda = 355\text{nm} \rightarrow S_{355} = S_{1064} / 0.85$<br>$\lambda = 10.6\mu\text{m} \rightarrow S_{10600} = S_{1064} / 3.5$<br>(10.6 $\mu\text{m}$ Reflectivity Warning) |      |
| 9  | Packing                  | Antistatic Bag With Label  | -    |

#### **4. Operating Instructions**

Operate within the specified laser damage thresholds.  
Keep clean and do not touch the active area.

#### **5. Electrical Connections and Filters**

It is recommended that both output channels are first filtered using 200 MHz low pass filters.  
The total output signal can be then extracted as the differential signal between the two signal cables.  
Output signal must be sampled with a minimum rate of 500 MSps.  
Peaks of typical signals range from few mVolts to few Volts.  
The calculated peaks areas are proportional to the pulse energy.  
Peak values are also proportional to pulse energy, but are dependent on energy power density, hence dependent on laser pulse duration.

#### **6. Warranty**

P/N BM-A-5W-14-TX-ENH is covered by one-year warranty from the date of shipment.  
Warranty applies against material and/or workmanship defects, provided that the instrument has been used under specified operating conditions; the warranty does not cover damages related to accident or misuse.  
During the warranty period Laser Point srl will repair or, at their option, replace any P/N BM-A-5W-14-TX-ENH or components that proves to be defective provided the parts are returned, shipping prepaid to Laser Point Customer Service or to another facility, authorized by Laser Point.  
Any attempt by an unauthorized person or entity to alter or repair the product voids its warranty.  
No other expressed warranty is given by Laser Point.  
Customers must fill in and mail the warranty card in order to validate the guarantee.  
In case of instrument malfunctioning or failure, contact Laser Point directly or its local distributor to obtain a Return Material Authorization Number (RMA).  
The material should be returned, transportation and insurance prepaid, to:

Laser Point srl  
Customer Service  
Via Burona, 51  
I-20055 Vimodrone  
Italy  
E-mail:  
[sales@LaserPoint.it](mailto:sales@LaserPoint.it)

Laser Point assumes no risk for the possible damage during shipping.

#### **7. Limited Liability**

For a free of charge warranty service, please provide a description of failure or encountered problems when filling the RMA request form.  
Laser Point will, at its option, repair or replace the defective product free of charge. However, if Laser Point determines that the failure was caused by misuse, alterations, accident or abnormal condition of operation or handling, Customer will be billed for the repair and the repaired product will be returned to Customer, transportation and insurance prepaid.

## **8. Compliance to RoHS Directive (RoHS 2015/863)**

The European (RoHS) Directive about Restriction of Hazardous Substances (RoHS 2015/863) aims to minimise the environmental impact of waste of electrical and electronic equipment by reducing the quantities of four heavy metals and two brominated flame retardants.

Laser Point P/N BM-A-5W-14-TX-ENH is a product compliant to RoHS European Directive.

## **9. European Union WEEE Directive (WEEE 2012/19/EU)**

The European Waste Electrical and Electronic Equipment (WEEE) Directive (WEEE 2012/19/EU) is represented by a crossed-out garbage container label (see below). The purpose of this directive is to avoid the disposal of WEEE as unsorted municipal waste and to facilitate its separate collection.



\*\*\*\*\*End of Document\*\*\*\*\*